


NOTEABLE™ Face Mask Aerosol Research at Baylor University November 18, 2020



Two experimental setups were tested with an aerosol generator and a detector: Setup A with a Noteable mask alone, and Setup B with a Noteable mask and choir robe positioned beneath.

For both setups, within the experimental parameters, no aerosol particles larger than 3 micrometers were detected in front of the mask, beneath the mask, or in front of the choir robe.

For more information, please contact:

Christie Sayes, PhD
Aerosol Scientist
Baylor University
Email: christie_sayes@baylor.edu